



Rapid Response Body Thermometer HY11P32 Application Note

Table of Contents

1	INTRODUCTION.....	2
2	THEORY DESCRIPTION.....	2
2.1	Sensor.....	2
2.2	CONTROL IC.....	3
2.3	MEASUREMENT THEORY.....	4
2.4	OUTPUT RATE.....	4
3	OPERATION DESCRIPTION.....	5
4	DESIGN SCHEME.....	6
4.1	HARDWARE DESCRIPTION.....	6
4.2	CIRCUIT DESCRIPTION.....	7
4.3	SOFTWARE DESCRIPTION.....	8
5	TECHNICAL SPECIFICATION.....	10
5.1	HARDWARE TECHNICAL SPECIFICATION.....	10
5.2	SOFTWARE TECHNICAL SPECIFICATION.....	10
6	CONCLUSION.....	10
7	REFERENCE.....	10
8	DEMON CODE.....	10

1 Introduction

One of the conventional methods to measure body temperature is by mercury thermometers. Heavy metal that leaked from broken mercury thermometers is harmful to the environment and human health; therefore electronic thermometers would be used to replace mercury thermometers gradually. The wide application of thermistor has created new direction for temperature measurement.

Based on the characteristic of different resistor value that Thermistor responses toward different temperature, we can easily realize temperature measurement. This article depicts how to apply HY11P32 of HYCON Technology to rapid response body temperature measurement.

2 Theory Description

2.1 Sensor

Thermistor is resistors that sense heat, which main function is to display resistor values that change with the environment. Thermistor is broadly divided into Positive Temperature Coefficient (PTC) and Negative Temperature Coefficient (NTC). NTC has high sensitivity and reasonable price, thus it was widely applied.

Users need to list a table before measurement, in which data can precisely reflect the conditions of Thermistor with changing temperature. For example, when temperature increases from 32.1°C to 32.2°C, the value of 503ET (NTC thermistor) drops from 36.763KΩ to 36.605KΩ. The temperature change is very limited, thus the process can be deemed as linearity. We table the value within human body temperature ranges (in 0.1°C temperature change, its resistor value is almost linear), using look-up-table method to gain its corresponding temperature.

NTC Thermistor 503ET equation is as follows:

Equation 2-1

$$R_{th}(T) = R_{25} \times \exp \left\{ B \times \left[\frac{1}{(T+273.13)} - \frac{1}{(25+273.13)} \right] \right\}$$

R_{th}(T): Thermistor change resistor value

B: Sensitivity of Thermistor

R₂₅: 25°C resistor value




Figure 1 Thermistor Resistor and Temperature Curve

2.2 Control IC

HY11P32 IC Features:

8-bit RISC-like simplified instruction set, 46 instructions included.

2.2V to 3.6V operation voltage, -40~85°C operation temperature range

Internal high precision RC oscillator, users can plan the most power efficient solution by 4 types of CPU operation clock source options.

Operation mode: 300μA@2MHz

Standby mode: 3μA@28KHz

Sleep Mode: 1μA

2K Word OTP (One Time Programmable) Type program memory, 128 Byte data memory

Brownout detector and Watch dog Timer, preventing CPU from crashing

18-bit fully differential input ΣΔADC (analog-to-digital converter)

Built-in PGA (Programmable Gain Amplifier), 10 signal amplifications, selecting from 1/4, 1/2, 1...to 128.

Built-in input zero point adjustment can add measurement range to different applications.

1.0V internal analog circuit common ground power source, equipped with Push-Pull drive power to provide driving voltage to sensors.

LVD low voltage detection function with 14 steps detection voltage configurations.

Analog power source, VDDA has 10mA regulated power output function.

4×12 LCD driver

1/4 Duty, 1/3 Bias

Built-in Charge Pump regulated circuit, providing 4 LCD bias.

8-bit Timer A

Build-In EPROM

3 Operation Description

1. Before countdown completed, enter into calibration mode by press PT1.0.
2. First WORD of EPROM is blank, entering into auto calibration mode.
3. Take 10 data and compare to the first data. If the difference falls in the range of 96H, start to judge. 10 data in successive falls in the range of 96H, then the last data would be deemed as calibration value to write EPROM otherwise, re-judge is needed.
4. When comparing, ADO is the judge of [(+) - (-)]/2 shift. Only when H=02H, M=01111XXXH will it enter into auto calibration judge mode to read value and store to EPROM.
5. If VPP does not connect to external 6V, then WAVPP message will show up when entering into calibration mode.
6. Using 37C reference resistor, 29938.5R to execute calibration, measurement will be restarted after calibration completed.
7. Temperature under 32C displays OVRL; higher than 43C display OVRH
8. Under measurement mode, pressing PT1.0 to start TIMERA, maximum value will be locked for around 4~5 seconds and display H xxxx, BZ ON
9. Then presses PT1.0 back to normal measurement mode
10. Under measurement mode, pressing PT1.1, displaying the maximum temperature value and display r xxxx
11. Then presses PT1.1, back to measurement mode
12. PT1.0/PT1.1 display maximum value and history maximum value does not function under OVRL/OVRH.
13. Under measurement mode, pressing PT1.2 to enter Sleep Mode, then press PT1.0 back to normal measurement mode

4 Design Scheme

4.1 Hardware Description




Figure 4 HY11P32 Simple Application Circuit

Main components:

1. MCU: HY11P32, main function is to measure electrical signal, control, operate and display
2. Sensor (TS): 503ET, to convert temperature and electrical signal
3. Reference resistance (RF): Reference resistor (precision resistor)

4.2 Circuit Description

Measuring temperature change by measuring the voltage change of Thermistor.

Reference resistor, R_F , is used to divide voltage. SD18 passes through R_F to gain a reference voltage. Accuracy of the measured resistor value relates to that of reference resistor. It is important to assure the accuracy of the reference resistor.

Calibration resistor, 29938.5Ω is the resistor value of Sensor 503ET under 37°C ambient temperature. When designing program, this value can be used for calibration.

BZ function is enabled after completion of locking max. value.

S0 has awakening IC function.

S1 has calibration, awakening IC and lock max. value function.

Press PT1.0 to enter calibration mode before the end of countdown when switched on.

Press PT1.0 to lock the max. value under measurement mode.

Press PT1.0 to awake the IC after entering into Sleep Mode.

S2 can display max. value in measurement history

S3 can enter into Sleep Mode

4.3 Software Description




Figure 5 Program Flow




Figure 6 Enter into Sleep Mode and Awaken Flow

5 Technical Specification

5.1 Hardware Technical Specification

Digital operation voltage: 2.2~3.6V

Analog operation voltage: 2.6~3.6V

Operation mode current (Measure Temperature 34.5°C): 404 μ A

Sleep mode current: 0.7 μ A

Temperature range: 32.1~42.5°C

Resolution: 0.01°C

Accuracy: 0.05°C

Data convert: 20Hz

5.2 Software Technical Specification

ADC configuration

ADC sampling rate: 250KHz

PGAGN: by pass PGA amplifier

ADGN: $\times 2$

OSR: 1024

ADC output rate: 250KHz/1024

Close SI \pm input buffer

Channel configuration

Senor signal

SI+: AI0

SI-: AI1

INH \rightarrow SI+, INL \rightarrow SI-

Voltage reference signal

VR+: AI2

VR-: VSS

Close VR \pm input Buffer

VREGN: $\times 1$

Power consumption control mode

Under measurement mode, if measurement temperature is within temperature range, press PT1.0 and enter to Standby mode.

Under measurement mode, if measurement temperature is within temperature range, press PT1.1 and enter to Standby mode.

Under measurement mode, press PT1.2 and enter to Sleep mode.

Under Sleep mode, press PT1.0 to awake IC

HY11P32 Application Note for Rapid Response Body Thermometer

6 Conclusion

Main features of using HY11P32 chip to realize rapid response temperature measurement solution:

1. Low current consumption
2. Fast output rate, 20 data/sec
3. Simple circuit and higher measurement precision
4. Reduce production line rework cost that caused by inaccurate temperature measurement
5. Gain better measurement data and stronger anti-interference performance in the same cost level.

Measurement precision and data is shown in below table:

IDEAL		ACTUAL		
TA(°C)	R(Ω)	TA(°C)	Rolling Counts	Deviation
42.5	24048	42.46	+/-1 count	-0.09%
42.0	24539	41.97	+/-1 count	-0.07%
41.0	25503	40.98	1 count	-0.05%
40.0	26517	39.99	1 count	-0.02%
39.0	27646	38.99	1 count	-0.03%
38.5	28186	38.49	1 count	-0.03%
38.0	28746	38.00	1 count	0.00%
37.5	29363	37.50	1 count	0.00%
37.0	29938	37.00	1 count	Calibration point
36.5	30573	36.50	1 count	0.00%
36.0	31194	36.01	1 count	0.03%
35.0	32557	35.01	1 count	0.03%
34.0	33924	34.02	1 count	0.06%
33.0	35386	33.02	1 count	0.06%
32.1	36763	32.13	1 count	0.09%

OSR:1024(244Hz); positive/negative network crosses to 12 data;
equivalent output rate is around 20Hz

Table 1 Measurement Deviation Table

7 Reference

[1] <http://www.hycontek.com/> HYCON Technology Corporation

8 Demon Code



HY11P32-Electronic Thermometer.rar